

Cylindrical torpedo												
												
Diamond Grit Size ●SC ●C ●F ●XF												
Dx Order No.	S217	S218	S316	S317								
size Ø 1/10 mm	010	012	012	014								
L mm	6.0	6.0	8.0	8.0								
Ref. No.	S877	S877	S878	S878								
ISO No.	288	288	289	289								

Flat end taper												
												
Diamond Grit Size ●SC ●C ●F ●XF												
Dx Order No.	S207	S208	S208B	S209B	S210	S211	S305	S306	S307	S310	S311	S313
size Ø 1/10 mm	009	010	012	016	010	014	012	014	014	012	014	018
L mm	3.0	4.0	4.0	4.0	7.0	7.0	8.0	8.0	8.0	10.0	10.0	10.0
Ref. No.	S845	S845	S845	S845	-	-	S847	S847	S847	S848	S848	S848
ISO No.	169	170	170	170	171	171	172	172	172	173	173	173

Round end taper												
												
Diamond Grit Size ●SC ●C ●F ●XF												
Dx Order No.	S251	S320A	S320	S321	S322	S323	S323B	S324B	S325			
size Ø 1/10 mm	012	010	012	014	016	018	018	014	016			
L mm	6.0	8.0	8.0	8.0	8.0	8.0	8.0	10.0	10.0			
Ref. No.	S855	S856	S856	S856	S856	S856	-	S850	S850			
ISO No.	197	198	198	198	198	198	198	199	199			

Conical torpedo												
												
Diamond Grit Size ●SC ●C ●F ●XF												
Dx Order No.	S221	S371	S372	S372B	S373							
size Ø 1/10 mm	014	012	014	016	016							
L mm	6.0	8.0	8.0	8.0	10.0							
Ref. No.	S877K	S878K	-	S878K	S879K							
ISO No.	297	298	298	298	299							

Bud												
												
Diamond Grit Size ●SC ●C ●F ●XF												
Dx Order No.	S300											
size Ø 1/10 mm	023											
L mm	5.0											
Ref. No.	S368											
ISO No.	257											

Pointed cone												
												
Diamond Grit Size ●SC ●C ●F ●XF												
Dx Order No.	S331	S332	S333	S335	S336	S337						
size Ø 1/10 mm	012	014	016	010	012	016						
L mm	8.0	8.0	8.0	10.0	10.0	10.0						
Ref. No.	S858	S858	S858	S859	S859	S859						
ISO No.	165	165	165	166	166	166						

Flame												
												
Diamond Grit Size ●SC ●C ●F ●XF												
Dx Order No.	S287	S340	S341	S342								
size Ø 1/10 mm	012	010	012	014								
L mm	5.0	8.0	8.0	8.0								
Ref. No.	S860	S862	S862	S862								
ISO No.	247	249	249	249								

Round edge wheel												
												
Diamond Grit Size ●SC ●C ●F ●XF												
Dx Order No.	S436	S516										
size Ø 1/10 mm	040	045										
L mm	1.3	1.5										
Ref. No.	S909	S909										
ISO No.	068	068										

FG Special

Guide pin												
												
Diamond Grit Size ●SC ●C ●F ●XF												
Dx Order No.	S36J	S38J										
size Ø 1/10 mm	018	018										
L mm	6.0	8.0										
Non-cutting mm	0.5	0.5										
Ref. No.	S80P	S81P										
ISO No.	-	-										

Round end taper with guide pin												
												
Diamond Grit Size ●SC ●C ●F ●XF												
Dx Order No.	S48J	S48J	S48K	S49J	S49K	S49L	S49MN					
size Ø 1/10 mm	016	018	021	016	021	023	026					
L mm	8.0	8.0	8.0	8.0	8.0	8.0	8.0					
Non-cutting mm	0.9	0.5	0.5	0.9	0.9	0.9	0.9					
Ref. No.	S86P	S86P	S86P	S98	S98	S98	S98					
ISO No.	508	508	508	508	508	508	508					

Flame with safe end												
												
Diamond Grit Size ●SC ●C ●F ●XF												
Dx Order No.	S56G	S60G										
size Ø 1/10 mm	012	012										
L mm	6.0	10.0										
Non-cutting mm	0.9	0.9										
Ref. No.	S61SE	S63SE										
ISO No.	255	256										

Depth maker												
												
Diamond Grit Size ●SC ●C ●F ●XF												
Dx Order No.	S79J	S79K	S80N									
size Ø 1/10 mm	016	021	027									
L mm	6.2	6.2	8.2									
Depth mm	0.3	0.5	0.8									
Ref. No.	S84	S84	S84									
ISO No.	552	552	552									