



DIAMOND BUR

DENTEX DENTAL INDUSTRIAL CORP.


8F, No. 8, Lane 348, Sec. 2, Chung Shang Rd.,
Chung-Ho Dist., New Taipei City 235, Taiwan.
E-mail: michelle.dentex@outlook.com
dentex@ms77.hinet.net
Tel: +886-2-22424333
Fax: +886-2-22424345
https://www.dentex.com.tw



Version: 2020.09.

RA Diamonds


Round



Diamond Grit Size
●SC ●C ●F ●XF

Dx Order No.	W101	W102	W103	W104	W105	W106	W360	W361
size Ø 1/10 mm	009	010	012	014	016	018	021	023
L mm	-	-	-	-	-	-	-	-
Ref. No.	W801	W801	W801	W801	W801	W801	W801	W801
ISO No.	001	001	001	001	001	001	001	001

Inverted cone



Diamond Grit Size
●SC ●C ●F ●XF

Dx Order No.	W121	W122	W123	W124	W125
size Ø 1/10 mm	009	010	012	014	016
L mm	1.1	1.1	1.3	1.6	1.6
Ref. No.	W805	W805	W805	W805	W805
ISO No.	010	010	010	010	010

Cylinder



Diamond Grit Size
●SC ●C ●F ●XF

Dx Order No.	W201B	W203	W204	W205	W206	W206B	W301	W302
size Ø 1/10 mm	009	010	012	014	012	014	012	014
L mm	3.0	4.0	4.0	4.0	6.0	6.0	8.0	8.0
Ref. No.	W835	W835	W835	W835	W836	W836	W837	W837
ISO No.	108	109	109	109	110	110	111	111

Flat end taper



Diamond Grit Size
●SC ●C ●F ●XF

Dx Order No.	W208	W208B	W209	W209B	W237
size Ø 1/10 mm	010	012	014	016	016
L mm	4.0	4.0	4.0	4.0	6.0
Ref. No.	W845	W845	W845	W845	W846
ISO No.	170	170	170	170	171

Pointed cone



Diamond Grit Size
●SC ●C ●F ●XF

Dx Order No.	W284H	W336	W337
size Ø 1/10 mm	014	012	016
L mm	6.0	10.0	10.0
Ref. No.	-	W859	W859
ISO No.	164	166	166

Bud



Diamond Grit Size
●SC ●C ●F ●XF

Dx Order No.	W296	W300
size Ø 1/10 mm	018	023
L mm	4.5	5.0
Ref. No.	W368	W368
ISO No.	257	257

Round end taper



Diamond Grit Size
●SC ●C ●F ●XF

Dx Order No.	W324	W324B	W325
size Ø 1/10 mm	012	014	016
L mm	10.0	10.0	10.0
Ref. No.	W850	W850	W850
ISO No.	199	199	199

HP Diamonds

Round



Diamond Grit Size
●SC ●C ●F ●XF

Dx Order No.	603	604	706
size Ø 1/10 mm	018	025	033
L mm	-	-	-
Ref. No.	H801	H801	H801
ISO No.	001	001	001

Round with collar



Diamond Grit Size
●SC ●C ●F ●XF

Dx Order No.	613	614	716
size Ø 1/10 mm	018	025	033
L mm	3.3	5.4	5.8
Ref. No.	H802	-	-
ISO No.	002	002	002

Inverted cone



Diamond Grit Size
●SC ●C ●F ●XF

Dx Order No.	634	636	637	803	903
size Ø 1/10 mm	016	021	025	040	060
L mm	1.9	1.9	2.4	1.4	2.5
Ref. No.	H805	H805	H805	H805	H805
ISO No.	010	010	010	014	014